

merident

MELECTIO



The light has been manufactured with a super-high luminosity LED. The intensity is from 15,000 to 30,000 Lux at 5 mm. These characteristics enable the light to support diagnosis during a therapeutic or surgical intervention.

- Super-high luminosity LED with assists in inspecting caries, calculus, fractures, and root canal orifices.
- Can reduce the need and expense of diagnostic X-rays.
- Immediate feedback saves time and improves patient care.
- Simple HI-LO-OFF one-button, battery-powered operation.



Visual imaging fracture, caries and root canals

Specifications

Battery type	2 x AAA alkaline batteries
Color temperature	White LED, 6500 K
Radiant intensity	Low mode:15,000 Lux above at 5 mm High mode:30,000 Lux above at 5 mm
Dimensions	L135 Ø 16 mm
Weight	Weight 59 g (including batteries)

Contents

- A** Melectio handpiece x1
- B** 3mm light guide rod x1
- C** AAA alkaline battery x2
- D** Instructions for use x1
- E** Disposable sleeve 20x
- F** Disposable Dental Mirror



Optional Accessories: Endo-Lux

3mm Light Guide Rod



Endo-Lux 5.5 Ø 1 with
10 pcs Endo Fiber



Endo-Lux Refill Pack



Disposable sleeves



Disposable Dental Mirror

